

# YMC CHIRAL

## Polysaccharide Derivatives Series

HPLC Column / Packing Material with Polysaccharide Derivatives Chiral Selector

- Applicable to various chiral compounds
- Excellent resolution/durability
- Extremely low initial cost on analysis and purification



# YMC CHIRAL

## Polysaccharide Derivatives Series

### Features



- Applicable to various chiral compounds
- Excellent resolution/durability
- Extremely low initial cost on analysis and purification
- High durability column that is suitable for SFC/SMB





Bulk packing material is also available.

### Specifications

#### Coated type

Column/Packing material	Particle size* <sup>1</sup> ( $\mu\text{m}$ )	Chiral selector	USP Classification
YMC CHIRAL Amylose-C	5 10 20	 Amylose tris(3,5-dimethylphenylcarbamate)	L51
YMC CHIRAL Cellulose-C		 Cellulose tris(3,5-dimethylphenylcarbamate)	L40
Usable mobile phase* <sup>3</sup>	<i>n</i> -hexane, <i>n</i> -heptane, ethanol, 2-propanol, acetonitrile, etc.		

#### Immobilized type\*<sup>2</sup>

Column/Packing material	Particle size* <sup>1</sup> ( $\mu\text{m}$ )	Chiral selector	USP Classification
YMC CHIRAL Amylose-SA	5 10 20	 Amylose tris(3,5-dimethylphenylcarbamate)	—
YMC CHIRAL Cellulose-SB		 Cellulose tris(3,5-dimethylphenylcarbamate)	—
Usable mobile phase* <sup>3</sup>	Normal phase	<i>n</i> -hexane, <i>n</i> -heptane, methanol, ethanol, 2-propanol, acetonitrile, ethyl acetate, tetrahydrofuran, chloroform, <i>t</i> -butyl methyl ether, etc.	
	Reversed phase	acetonitrile, methanol, ethanol, 2-propanol, tetrahydrofuran, water, aqueous buffer, etc.	

\*<sup>1</sup> Please inquire us other than those listed above.\*<sup>2</sup> Other immobilized type products will be launched.\*<sup>3</sup> Please also refer to the instruction manual available on our website.




## Useful for Chiral Separation of Wide Range of Compounds

Compound	Mobile phase	Separation factor ( $\alpha$ )							
		Coated type				Immobilized type			
		Amylose-C	Competitor's product	Cellulose-C	Competitor's product	Amylose-SA	Competitor's product	Cellulose-SB	Competitor's product
<i>trans</i> -Stilbene oxide	Hex/IPA (90/10)	2.9	3.0	2.3	2.2	2.7	2.8	1.6	1.9
Benzoin	Hex/IPA (90/10)	1.3	1.3	1.6	1.6	1.2	1.2	1.4	1.4
<i>N</i> -CBZ-DL-Alanine	Hex/IPA/TFA (80/20/0.1)	2.0	2.2	3.0	2.9	1.7	1.7	1.7	1.8
Ibuprofen	Hex/IPA/TFA (99/1/0.1)	1.1	1.1	1.3	1.2	1.1	1.1	1.1	1.1
Propranolol	Hex/IPA/DEA (80/20/0.1)	×	×	2.0	1.8	×	×	1.6	1.4
Verapamil	Hex/IPA/DEA (90/10/0.1)	1.3	1.3	×	×	1.2	1.2	×	×


Hex: *n*-hexane, IPA: 2-propanol, TFA: trifluoroacetic acid, DEA: diethylamine, ×: Not separated

YMC Chiral Polysaccharide Derivatives Series provide results comparable to other polysaccharide columns.

## Excellent Peak Shape




Propranolol



Column : 5  $\mu$ m, 250 X 4.6 mm I.D.  
Eluent : *n*-hexane/2-propanol/diethylamine (80/20/0.1)  
Flow rate : 1.0 mL/min  
Temperature : 25°C  
Detection : 230 nm  
Injection : 10  $\mu$ L (0.1 mg/mL)

YMC CHIRAL Polysaccharide Derivatives Series provide good peak shapes on ionic and metal coordination compounds.


## Low Column Bleeding



Gradient test  
Column : 5  $\mu$ m, 250 X 4.6 mm I.D.  
Eluent : A) *n*-hexane  
          B) ethanol  
          2-80%B (0-30 min)  
Flow rate : 1.0 mL/min  
Temperature : 25°C  
Detection : UV at 230 nm

YMC CHIRAL Amylose-SA and Cellulose-SB show remarkably reduced background signal under the typical gradient condition. This low column bleeding of those columns provides high sensitivity on LC/MS analysis due to the very low ion suppression as well as stable baseline. YMC CHIRAL Amylose-SA and Cellulose-SB offer excellent robustness on gradient analysis and highly sensitive analysis on LC/MS.

## Extended Packing Durability





### Sequential gradient test

Column : 5  $\mu$ m, 250 X 4.6 mmI.D.  
 Eluent : A) *n*-hexane, B) ethanol  
 0-100%B (0-15 min)  
 Flow rate : 3.0 mL/min  
 Pressure : 10-30 MPa/run  
 Temperature : 37°C


### Column performance test

Column : 5  $\mu$ m, 250 X 4.6 mmI.D.  
 Eluent : *n*-hexane/ethanol (90/10)  
 Flow rate : 1.0 mL/min  
 Temperature : 37°C  
 Detection : UV at 230 nm  
 Sample : *trans*-Stilbene oxide



YMC CHIRAL Polysaccharide Derivatives Series have outstanding packed bed stability provided by the use of high-strength super wide pore silica and innovative packing technology. The column efficiency and selectivity are maintained even after the sequential gradient tests at a high flow rate (three times higher than normal flow rate) and a high pressure (rapid pressure change). YMC CHIRAL Polysaccharide Derivatives Series are suitable for the case where one intends to shorten an analysis time, (re-)equilibration time and/or column cleaning time by increasing the flow rate. This feature is also effective when highly viscous solvent is needed to be used as a mobile phase on immobilized type columns.

## Wide Usable pH Range





### Continuous flow of acid/alkaline solution

Column : YMC CHIRAL Cellulose-SB  
 5  $\mu$ m, 50 X 4.6 mmI.D.  
 Eluent : Buffer/methanol (90/10)  
 Flow rate : 1.0 mL/min  
 [Acidic condition]  
 Buffer : 0.1% H<sub>3</sub>PO<sub>4</sub> (pH 2)  
 Temperature : 40°C  
 [Basic condition]  
 Buffer : 20 mM NH<sub>4</sub>HCO<sub>3</sub>-DEA (pH 9)  
 Temperature : 25°C

### Column performance test


Column : YMC CHIRAL Cellulose-SB  
 5  $\mu$ m, 50 X 4.6 mmI.D.  
 Eluent : acetonitrile/water (30/70)  
 Flow rate : 1.0 mL/min  
 Temperature : 25°C  
 Detection : UV at 254 nm  
 Sample : Benzoin



YMC CHIRAL Amylose-SA and Cellulose-SB have excellent chemical stability and can be used across a wide pH range. Addition of acid/amine in a mobile phase or use of buffer solution is an effective tool for improving peak shape and/or resolution.

## High Solvent Versatility (Immobilized Type)

### High Solvent Resistance



### Retention rate of initial column performance

(after flushing with 1,000CV of each solvent at 40°C)

\*CV=Column Volume

	Amylose-SA		Cellulose-SB	
	$\alpha$	$k'(2)$	$\alpha$	$k'(2)$
Ethyl acetate	100.3%	101.2%	100.0%	99.1%
Tetrahydrofuran	100.0%	100.0%	99.3%	98.0%
Dichloromethane	100.3%	100.6%	101.3%	99.6%

Column : 5  $\mu$ m, 50 X 4.6 mm I.D.  
Eluent : *n*-hexane/2-propanol (95/5)  
Flow rate : 1.0 mL/min  
Temperature : 25°C  
Sample : Benzoin



On YMC CHIRAL Amylose-SA and Cellulose-SB, the change in column performance after flushing with each solvent was less than 2%. They have high resistance to various solvents.

### Method Screening

Solvent Range on Coated Type

Solvent Range on Immobilized Type

*n*-hexane/ethyl acetate or MTBE/R-OH is also effective.





2-Phenylcyclohexanone

Column : YMC CHIRAL Cellulose-SB  
5  $\mu$ m, 150 X 3.0 mm I.D.  
Flow rate : 0.425 mL/min  
Detection : UV at 220 nm  
Temperature : 25°C


The chiral method screening of 2-Phenylcyclohexanone on YMC CHIRAL Cellulose-SB is shown above. A mobile phase containing MTBE gave good separation. On YMC CHIRAL Amylose-SA and Cellulose-SB with high solvent versatility, chromatographers can freely choose the most suitable mobile phase by considering the solubility, resolution and loadability of target compound based on the purpose of separation (e.g. analytical or preparative).

## Effective for Preparative Separation of Enantiomers

### Suitable for High-Loading Preparation



Increasing  
loading  
X 20



Competitor's product

YMC CHIRAL  
Cellulose-C



Competitor's product

YMC CHIRAL  
Cellulose-C

On both analytical and preparative separation, separation behavior of YMC CHIRAL Polysaccharide Derivatives Series is equivalent to that of the competitor's product.

Column : 5  $\mu$ m, 250 X 4.6 mm.I.D.  
Eluent : *n*-hexane/ethanol (90/10)  
Flow rate : 1.0 mL/min

### High Purity Purification Utilizing Recycling Preparation



	Single Column	Recycle
Enantiomeric purity (%ee)		
Fr.1 (R)-(-)-Ibuprofen	95.0	98.4
Fr.2 (S)-(+)-Ibuprofen	96.8	99.2
Yield (%)	84	95
Solvent consumption (mL solvent/g product)	9,523	1,276


Ibuprofen enantiomers were purified by utilizing recycling mode of multi preparative HPLC system, LC-Forte/R. Recycling chromatography is effective when method optimization of chiral isolation is difficult on single column separation. Recycling method offers high purity and high recovery purification. Furthermore, no solvent is consumed during recycling mode. It greatly contributes reduction of solvent consumption on purification.



LC-Forte/R

## Applications


### Valsartan (USP method)



Column : YMC CHIRAL Cellulose-C 5  $\mu$ m, 250 X 4.6 mmI.D.  
 Eluent : *n*-hexane/2-propanol/TFA (85/15/0.1)  
 Flow rate : 0.8 mL/min  
 Temperature : 25°C  
 Detection : UV at 230 nm

E130620B


### Sertraline (EP method)



Column : YMC CHIRAL Amylose-C 5  $\mu$ m, 250 X 4.6 mmI.D.  
 Eluent : mixture\*/*n*-hexane (70/30)  
 \**n*-hexane/2-propanol/diethylamine (975/25/1)  
 Flow rate : 0.4 mL/min  
 Temperature : 25°C  
 Detection : UV at 275 nm

E131205C


### Duloxetine (USP method)



Column : YMC CHIRAL Cellulose-C 5  $\mu$ m, 250 X 4.6 mmI.D.  
 Eluent : *n*-hexane/2-propanol/diethylamine (83/17/0.2)  
 Flow rate : 1.0 mL/min  
 Temperature : 40°C  
 Detection : UV at 230 nm

F130930A


### Verapamil



Column : YMC CHIRAL Amylose-C 5  $\mu$ m, 250 X 4.6 mmI.D.  
 Eluent : *n*-hexane/2-propanol/diethylamine (90/10/0.1)  
 Flow rate : 1.0 mL/min  
 Temperature : 25°C  
 Detection : UV at 254 nm

V130905D


### Donepezil



Column : YMC CHIRAL Amylose-SA 5  $\mu$ m, 250 X 4.6 mmI.D.  
 Eluent : *n*-hexane/2-propanol/diethylamine (80/20/0.1)  
 Flow rate : 1.0 mL/min  
 Temperature : 25°C  
 Detection : UV at 230 nm

V140422B2

### Astaxanthin



Courtesy of Fuji Chemical Industry Co., Ltd.

Column : YMC CHIRAL Cellulose-SB 5  $\mu$ m, 250 X 4.6 mmI.D.  
 Eluent : *n*-hexane/THF (85/15)  
 Flow rate : 0.5 mL/min  
 Temperature : 25°C  
 Detection : UV at 476 nm

E140417B

## Ordering Information

### Columns

Particle size ( $\mu\text{m}$ )	Column size inner diameter X length (mm)	Product number			
		Coated type		Immobilized type	
		Amylose-C	Cellulose-C	Amylose-SA	Cellulose-SB
5	4.6 X 150	KAN99S05-1546WT	KCN99S05-1546WT	KSA99S05-1546WT	KSB99S05-1546WT
	4.6 X 250	KAN99S05-2546WT	KCN99S05-2546WT	KSA99S05-2546WT	KSB99S05-2546WT
	10 X 250	KAN99S05-2510WT	KCN99S05-2510WT	KSA99S05-2510WT	KSB99S05-2510WT
	20 X 250	KAN99S05-2520WX	KCN99S05-2520WX	KSA99S05-2520WX	KSB99S05-2520WX

Please inquire us other than those listed above.


### Packing Materials

Particle size ( $\mu\text{m}$ )	Product number			
	Coated type		Immobilized type	
	Amylose-C	Cellulose-C	Amylose-SA	Cellulose-SB
5	KAN99S05	KCN99S05	KSA99S05	KSB99S05
10	KAN99S11	KCN99S11	KSA99S11	KSB99S11
20	KAN99S21	KCN99S21	KSA99S21	KSB99S21

Please inquire us other than those listed above.

## Contract Purification of Chiral Compounds is Now Available

YMC CO., LTD. Chiral Technologies Laboratory was established in November, 2013 and YMC has started contract purification of chiral compounds. Following the steps below, YMC is confident that highly-purified chiral compounds are delivered with high credibility and speed. YMC's contract purification service enables you to accelerate the realization of your research and development.



### Worldwide Availability

**YMC America, Inc.**  
www.ymcamerica.com

**YMC Europe GmbH**  
www.ymc-europe.com

**YMC India Pvt.Ltd.**  
www.ymcindia.com

**YMC Korea Co., Ltd.**  
www.ymckorea.com

**YMC Shanghai Rep. Office**  
www.ymchina.com

**YMC Singapore Tradelinks Pte. Ltd.**  
www.ymc.sg

**YMC Taiwan Co., Ltd.**  
www.ymctaiwan.com



### YMC CO., LTD.

YMC Karasuma-Gojo Bldg., 284 Daigo-cho  
Karasuma Nishiiru Gojo-dori, Shimogyo-ku  
Kyoto, 600-8106, Japan  
TEL: +81-75-342-4515 FAX: +81-75-342-4550  
www.ymc.co.jp sales@ymc.co.jp

Distributor