


## High resolution analysis of Oligonucleotides on reversed phase chromatography

F121018AE


### Oligonucleotides analysis on reversed-phase ion-pair chromatography (RP-IPC)

#### Comparison of retention and separation under various mobile phase conditions




### Applicability to LC/MS analysis

#### Impact of concentration and types of ion-pairing reagent on resolution and signal intensity



### LC/MS analysis of miRNA




Courtesy of M. Yamada, SHIMADZU CORPORATION

# High temperature analysis of oligonucleotides with YMC-Triart C18


## Effect of mobile phase and column temperature on separation of siRNA duplex

Crude synthetic siRNA duplex (19 bp) : 5'-CGU ACG CGG AAU ACU UCG AdTdT-3'  
3'-dTdTGCA UGC GCC UUA UGA AGC U-5'



- Separation of siRNA duplex under different mobile phase conditions at various temperatures with YMC-Triart C18 is shown.
- Under both condition A and condition B, peak shape and resolution between immediate peaks is improved by increasing the column temperature.
- Due to the improvement of dispersion and distribution velocity when increasing column temperature, bio-macro molecules such as RNA and DNA generally exhibit sharper peak shape and improved resolution.
- Under condition B at 40 °C or higher temperature, two peaks of single-stranded RNA that is generated by denaturation of siRNA duplex are observed. This HPLC technique that is utilizing high temperature to generate single-stranded RNA is called "Denaturing HPLC", and widely used in the field of gene mutation analysis.
- As shown above, denaturation of duplex DNA or RNA is also influenced by ionic strength (type and concentration), pH and polarity as well as temperature. Those analysis conditions (temperature and mobile phase) are recommended to be optimized depending on characteristics of target analyte and purpose of analysis.

## Durability at pH 6.0 (DBAA buffer) and 65°C



<b>Test condition</b>	Column : 1.9 $\mu$ m or 2.0 $\mu$ m, 12 nm, 50 X 2.0 mm I.D.
	Eluent : A) 10 mM di-n-butylamine-acetic acid (pH 6.0)
	B) methanol
	30-50%B (0-20 min)
	Flow rate : 0.4 mL/min
	Detection : UV at 269 nm
	Temperature : <b>65°C</b>
	Sample : Oligodeoxythymidylic acid, [d(T) <sub>2-20</sub> ]
	Injection : 1 $\mu$ L (5 nmol/mL)
	System : Agilent 1290

- Combination of neutral buffer containing amino ion-pairing reagent and high temperature is useful for high-throughput analysis of oligonucleotides or denaturing HPLC. However, conventional silica-based reverse-phase column can hardly be used with such condition due to the poor durability.
- YMC-Triart C18 using inorganic/organic hybrid silica with thorough surface modification offers excellent durability at elevated temperature and pH. YMC-Triart C18 is ideal for oligonucleotides analysis.