



Enantiomeric Purification Utilizing Supercritical Fluid Chromatography (SFC)

Purification of Flavanone
Using Alcyon SFC CSP Amylose-C

YMC CO., LTD.

Features of Alcyon SFC column

- Available in chiral and achiral stationary phases
- Faster separation with high resolution
- Excellent durability
- Great reduction of solvent consumption



SFC preparative separation of Flavanone


Conditions

HPLC

Column : CHIRAL ART Amylose-C
5 μ m, 250X20 mmI.D.
Eluent : *n*-hexane/ethanol (90/10)
Flow rate : 20 mL/min
Detection : UV at 220 nm
Injection : 3 mL (20 mg/mL)
Temperature : ambient

SFC

Column : Alcyon SFC CSP Amylose-C
5 μ m, 250X20 mmI.D.
Eluent : CO₂/ethanol (80/20)
Flow rate : 60 mL/min
Detection : UV at 280 nm
Injection : 1.5 mL (20 mg/mL)
Temperature : 30
Back Pressure : 15 MPa




Flavanone

SFC preparative separation of Flavanone




Preparative chromatograms

HPLC



HPLC	Fr.1	Fr.2
Enantiomeric purity	100 % ee	99.7 % ee
Yield	95.7%	93.7%

SFC



SFC	Fr.1	Fr.2
Enantiomeric purity	100 % ee	99.8 % ee
Yield	94.5%	95.6%

Cost effective purification using SFC

	HPLC		SFC	
Column : 250 X 20mmI.D.	Fr.1	Fr.2	Fr.1	Fr.2
Enantiomeric purity (%ee)	>99.9	99.7	>99.9	99.8
Yield (%)	95.7	93.7	94.5	95.6
Productivity* ¹ (mg product / hr)	172	169	340	344
Fractionated liquid volume (L solvent / g product)	1.15	2.88	0.39	0.57
Solvent consumption (L solvent / g product)	about 7		about 2	
Cost factor* ²	1		0.417	

*1 Calculated based on injections of every 2.5 minutes on SFC, and of every 9.0 minutes on HPLC.

*2 Calculated based on costs of solvents/CO₂ and waste disposal. Cost on SFC is calculated when the cost on HPLC is fixed as 1.

Effective of SFC preparative method compared to HPLC preparative method

- Productivity of SFC is twice as high as that of HPLC.
- Fraction volume on SFC is reduced by 33%, and as a result, solvent is easily removed at post-purification process.
- Solvent consumption on SFC is reduced by 25% SFC.
- Total purification cost on SFC is less than half of that on HPLC purification.