


## Efficient method development for chiral separation by using CHIRAL ART columns

U190520AE


Liquid chromatography using polysaccharide-based chiral separation columns is a popular method for the separation of enantiomers. However, it has been very hard to find the correct chiral separation column and mobile phase condition for achieving an acceptable chiral separation.

CHIRAL ART columns (including four immobilized types and two coated types) exhibit excellent chiral separation performance for a wide range of racemic compounds. In this report, we introduce the efficient method development for the separation of enantiomers via a high throughput method scouting using the CHIRAL ART columns.

### Workflow of method development for chiral separation




### Method scouting strategy for chiral separation



## Example of method scouting by using CHIRAL ART columns

The method scouting for the chiral separation of phenoxybenzamine was carried out. As shown in below chromatograms, the baseline separation of phenoxybenzamine was achieved using Cellulose-SJ column. For the further method development, the combination of Cellulose-SJ and *n*-hexane/ethanol containing 0.1% DEA was chosen according to resolution and peak shape.



## Optimization of chromatographic conditions


On the basis of the scouting results, the chromatographic conditions were optimized depended on the purpose of the chromatography experiment. To operate in an isocratic elution mode, it is appropriate to decrease the proportion of the polar organic solvent component in mobile phases by 10-15%pt from the proportion at the first eluting peak under gradient conditions. For a rapid analysis, short columns with 3  $\mu$ m particles are recommended. For a preparative purpose, the method optimization and loading studies by using 250 X 4.6 mmI.D. columns with 5  $\mu$ m particles make it easy to scale up to preparative columns.

### Optimized conditions


- Eluent (gradient  $\Rightarrow$  isocratic)
- Particle size, Column size
- Flow rate
- Injection volume, etc.

### Chosen chromatogram via the method scouting

Column : CHIRAL ART Cellulose-SJ  
 3  $\mu$ m, 75 X 3.0 mmI.D.  
 Eluent : A) *n*-hexane/DEA (100/0.1)  
 B) ethanol/DEA (100/0.1)  
 5-45%B (0-6 min)



Optimization  $\rightarrow$



Column : CHIRAL ART Cellulose-SJ  
 5  $\mu$ m, 250 X 4.6 mmI.D.  
 Eluent : *n*-hexane/ethanol/DEA (95/5/0.1)  
 Flow rate : 1.0 mL/min  
 Temperature : 25 °C  
 Detection : UV at 270 nm  
 Injection : 5  $\mu$ L (1.0 mg/mL)