

High resolution analysis of monoclonal antibodies and fragments using YMC-Pack Diol

R170831AE


YMC-Pack Diol is a silica gel based size exclusion chromatography (SEC) phase which comes in range of particle sizes and pore sizes for the optimum separation of different biomolecules. Columns packed with 2 μ m particles are effective for rapid analysis and high resolution of aggregates and fragments of biopharmaceuticals.

Analysis of monoclonal antibody and its aggregates

● Influence of flow rate on antibody analysis

YMC-Pack Diol-300, 2 μ m, 300 X 4.6 mmI.D.

Flow rate




Flow rate (mL/min)	N (2)	P (MPa)
0.8	7,300	28.3
0.6	8,900	20.9
0.4	9,400	14.4
0.2	15,200	7.2

Column : YMC-Pack Diol-300, 2 μ m, 300 X 4.6 mmI.D.
 Eluent : 0.1 M KH_2PO_4 - K_2HPO_4 (pH 7.0) containing 0.2 M NaCl
 Detection : UV at 280 nm
 Temperature : ambient
 Sample : Humanized monoclonal IgG1

- Diol-300 2 μ m column was used for monoclonal antibody analysis. Resolution between aggregates and monomer is improved as the flow rate is decreased.
- Higher flow rates are suitable for increased sample throughput, especially when the resolution is sufficient.

● Increasing throughput by using the 2 μ m column




Particle size (μ m)	Column size length X I.D. (mm)	N (3)	Rs (1,2)	Rs (2,3)
5	300 X 4.6	8,500	0.88	2.67
2	300 X 4.6	16,200	1.17	4.15
2	150 X 4.6	8,700	0.85	2.75



Column : YMC-Pack Diol-300, 2 μ m,
 Eluent : 0.1 M KH_2PO_4 - K_2HPO_4 (pH 7.0) containing 0.2 M NaCl
 Flow rate : 0.2 mL/min
 Detection : UV at 280 nm
 Temperature : ambient
 Sample : Humanized monoclonal IgG1

- The 2 μ m, 150 mm length column offers the same resolution as the 5 μ m, 300 mm length column. This means the analysis time can be reduced by half by changing the particle size from 5 μ m to 2 μ m.

Analysis of monoclonal antibody fragments by using YMC-Pack Diol-200 column




Method of papain digestion



Column	: YMC-Pack Diol-200, 2 μm 300 X 4.6 mm I.D.
Eluent	: 0.1 M KH ₂ PO ₄ -K ₂ HPO ₄ (pH 7.0) containing 0.2 M NaCl
Flow rate	: 0.2 mL/min
Detection	: UV at 280 nm
Temperature	: ambient

- A monoclonal antibody digested by papain, a proteolytic enzyme, as well as the intact monoclonal antibody were analyzed using a YMC-Pack Diol-200 column.
- By using a high resolution 2 μm column, intact antibody and fragments were well separated.

Lot-to-lot reproducibility of YMC-Pack Diol 2 μm



Column	: YMC-Pack Diol-300, 2 μm 150 X 4.6 mm I.D.
Eluent	: 0.1 M KH ₂ PO ₄ -K ₂ HPO ₄ (pH 7.0) containing 0.2 M NaCl
Flow rate	: 0.2 mL/min
Detection	: UV at 280 nm
Temperature	: 25°C
Sample	: Humanized monoclonal antibody

- YMC-Pack Diol 2 μm columns offer high lot-to-lot separation reproducibility. This feature greatly contributes to characterization of antibodies in areas of research as well as quality control of biopharmaceuticals.