



May 09, 2024

To Valued Customers

YMC CO., LTD.
Quality Assurance Department

Notice Regarding BioPro IEX Columns

We are planning to upgrade the quality control chromatography system we use for shipping control of BioPro IEX columns as it is obsolete. We will change the COA format and serial number to reflect the change. In addition, the airtight stopper change to PEEK, the same as the column material.

Please kindly take a look at this notification. If you should have any questions or concerns in this change, please do not hesitate to contact us. We thank you for your understanding.

1. Product to be changed
BioPro IEX columns
2. Date of change
Change sequentially from shipments in and after June 2024
3. Contact information
Website Inquiry Form: <https://www.ymc.co.jp/en/inquiry/form/>



Appendix

Details of Changes to COA format and Serial number

	Current	New
COA format		
Serial number	BP○○○○○○○ e.g. BP123456	○○○□□○○○○○ e.g. 123AB12345
	○ : Number, □ : Alphabet	

Details of change to airtight stopper

	Current	New
Material	POM	PEEK
Color	Black	Brown